


Tested to ensure  
non-pyrogenicity (LAL), HSSA  
and sperm compatibility



## BENEFITS

### Shapeability

Shapeable IUI catheters are manufactured with superior memory enabling the doctor to form the tip of the catheter to align with the cervical canal

Clear flexible & shapeable IUI catheters allow for full visibility of the washed specimen as it is injected into the uterus. The clear shaft also reduces the possibility of uterine spasms and retroflow by removing any visible air bubbles in catheter.

### Directional Injection

Side hole near the tip is aligned with the printed depth marks and a flat side of the connector. The alignment allows the catheter to point towards either the left or right fallopian tube ostium for directional injection.

Download IFU



## TECHNICAL DATA

French Size . . . . . 5Fr

Dead Space Volume. . . . . see ordering information (ml)

Luer Type . . . . . Female

Box Qty. . . . . 25 ea.



Not made with natural rubber latex.



**ORDERING INFORMATION**

Deadspace volumes listed below (.XXml)

- TMI1192..... 15cm, 5Fr IUI Catheter (.25ml)
- TMI1193..... 15cm, 5Fr Shapeable IUI Catheter (.15ml)
- TMI1194..... 18cm, 5Fr IUI Catheter (.28ml)
- TMI1195..... 18cm, 5Fr Shapeable IUI Catheter (.20ml)
- TMI1295..... 18cm, 5Fr Clear Shapeable IUI Catheter (.20ml)
- TMI1196..... 23cm, 5Fr IUI Catheter (.32ml)
- TMI1197..... 23cm, 5Fr Shapeable IUI Catheter (.25ml)
- TMI1297..... 23cm, 5Fr Clear Shapeable IUI Catheter (.25ml)



**Thomas Medical is an industry leading manufacturer of single-use women's health products.**

[www.thomasmedical.com](http://www.thomasmedical.com) | (800) 556-0349 | [sales@thomasmedical.com](mailto:sales@thomasmedical.com)  
6102 Victory Way Indianapolis, IN 46278

**CE** **Made in the USA**  
ISO 13485 and FDA Registered